

CARE Myocardium 失败取证 Deep Research 证据包

CARE Forensic Research Controller

20260730 本地证据冻结版

CARE Myocardium 失败取证 Deep Research 证据包

本 PDF 使用 Pandoc + XeLaTeX final-standard 路线生成, 拉丁字体为 TeX Gyre Termes, 中文字体来自 /users/a/e/aereinh/render_resources/chinese_math_pdf 中的 NotoSerifSC/NotoSansSC。它不是新模型蓝图, 不包含 validation upload, 也不声明 hosted 指标。

一页执行摘要

V2 的实际结论是: 历史 CARE 路线长期未稳定超过 nnU-Net, 主要不是某一个概念天然错误, 而是强基线继承、decoder 完整性、final-mask 组件进入路径、病例级 help/harm 选择、标签/评价语义和训练/recipe 绑定没有同时闭合。V2 已补齐 G1-G10 的终态证据; 其中缺 exact asset 的项目按 BLOCKED_BY_MISSING_BOUND_ASSET 写入, 不再把缺失证据伪装成负结果。




Figure 1: 证据等级计数

目录式章节索引

章节	主题
1	为什么现在必须做失败取证
2	CARE 数据、中心、模态和标签真值
3	官方与内部指标语义
4	当前评价代码中的已确认问题
5	nnU-Net 强基线到底强在哪里
6	SRR v2-v3 的设计意图与落地差距
7	Batch 0-7 历史证据
8	MMRD 的设计、实现和失败
9	Cascade/DG 的设计、实现和失败
10	ARC 的设计、实现和失败
11	PRISM W1-W3 的完整复盘
12	MoSAIC clean、full-data 和 hosted recipe
13	所有模型统一病例级比较
14	困难子组
15	case-wise help/harm
16	失败病例视觉图册
17	错误重合和模型互补上限
18	selector feasibility
19	冻结特征可分性 probe
20	decoder-reset 诊断对照
21	多序列错位是否为主因
22	scar 的真实瓶颈
23	pure edema 的真实瓶颈
24	Cine 的真实瓶颈
25	为什么过去多次充分设计仍然失败
26	根因排序与证据图
27	当前能下的结论
28	当前不能下的结论
29	外部 Deep Research 必须回答的问题
30	下一轮决策树

1. 为什么现在必须做失败取证

过去几轮路线没有稳定超过 nnU-Net, 不能直接归结为“模型不够复杂”。更可靠的取证路径是把设计承诺、实现连线、训练预算、checkpoint 选择、评价语义、预测缓存和 hosted recipe 分开冻结。

2. CARE 数据、中心、模态和标签真值

本节回答数据层面是否存在足够明确的标签和模态条件。关键边界是 official scar、official pure edema 和 internal edema-zone 必须分开。




Figure 2: 中心病例数




Figure 3: 病灶体积分布

cases	scar_positive	pure_edema_positive	t2_present
220	212	80	220

3. 官方与内部指标语义

第 3 页不使用宽表。下面只列三列：对象、内部标签、允许声明范围，避免右侧列截断。

object	internal_labels	allowed_claim_scope
scar	5	official
pure_edema	4	T2-present official edema
edema_zone	4 5	internal only

reference evaluator 的 known-bad fixtures 覆盖 remote FP、spacing HD95、empty case、lesion recall 和 label 4/5 语义。V2 对可绑定预测执行统一病例级重聚合；缺 exact asset 的旧模型不写成科学负结果。

4. 当前评价代码中的已确认问题

当前可确认的是评价风险，而不是所有历史结论已经被推翻。需要重算的对象包括 remote FP、HD95 physical spacing、empty-GT population mean，以及 pure edema 与 edema-zone 的混写。

5. nnU-Net 强基线到底强在哪里

nnU-Net 作为强基线的意义在于完整 decoder、稳定训练 recipe、成熟数据增强和直接 final mask 输出。当前包没有使用 foreground mean 掩盖 scar/pure edema。

6-10. SRR、Batch、MMRD、Cascade/DG、ARC 的历史证据

这些路线的历史证据等级不能混用。V2 将 Batch0-7、MMRD、Cascade、ARC、DG/DR/DPR 与 PRISM 分别绑定 source、checkpoint、prediction、metric 和 controller packet；缺 exact replay 资产的项目保持阻塞状态。

model_id	checkpoint_files_bound	prediction_files_bound	metric_files_bound	terminal_status
BATCH0_3_SRR_V2_3	ANCHOR_CONTROL	132	1	COMPLETED_WITH_VALID_
BATCH7_BR2_SIP	21	0	0	BLOCKED_BY_MISSING_BO
MMRD_BATCH9	9	0	0	BLOCKED_BY_MISSING_BO
SRR_CASCADE_RESCUE		0	2	COMPLETED_WITH_VALID_
CARE_ARC	8	0	1	COMPLETED_WITH_VALID_
CARE_DG_DR_DPR	129	132	2	COMPLETED_WITH_VALID_
CARE_PRISM_V2	13	455	1	COMPLETED_WITH_VALID_

- **NNUNET**
 - result_evidence_grade: A_VERIFIED_FAIR_FINAL_MASK
 - current_scientific_conclusion: 强基线；需继续绑定五折和同口径病例级指标。
- **SRR_V2**
 - result_evidence_grade: E_STALE_OR_INCONSISTENT
 - current_scientific_conclusion: 历史证据需绑定代码、checkpoint、split 和预测。
- **SRR_V25**
 - result_evidence_grade: E_STALE_OR_INCONSISTENT
 - current_scientific_conclusion: 历史证据需绑定代码、checkpoint、split 和预测。
- **SRR_V3**
 - result_evidence_grade: E_STALE_OR_INCONSISTENT
 - current_scientific_conclusion: 历史证据需绑定代码、checkpoint、split 和预测。
- **BATCH0**
 - result_evidence_grade: E_STALE_OR_INCONSISTENT
 - current_scientific_conclusion: 历史证据需绑定代码、checkpoint、split 和预测。
- **BATCH1**
 - result_evidence_grade: E_STALE_OR_INCONSISTENT
 - current_scientific_conclusion: 历史证据需绑定代码、checkpoint、split 和预测。
- **BATCH2**
 - result_evidence_grade: E_STALE_OR_INCONSISTENT
 - current_scientific_conclusion: 历史证据需绑定代码、checkpoint、split 和预测。
- **BATCH3**
 - result_evidence_grade: E_STALE_OR_INCONSISTENT
 - current_scientific_conclusion: 历史证据需绑定代码、checkpoint、split 和预测。
- **BATCH4**
 - result_evidence_grade: E_STALE_OR_INCONSISTENT
 - current_scientific_conclusion: 历史证据需绑定代码、checkpoint、split 和预测。
- **BATCH5**
 - result_evidence_grade: E_STALE_OR_INCONSISTENT
 - current_scientific_conclusion: 历史证据需绑定代码、checkpoint、split 和预测。

11. PRISM W1-W3 的完整复盘

PRISM 不能只看是否有强 encoder。V2 已完成 13 checkpoint replay 和 D0-D3 decoder-reset 诊断。最关键的负证据是：完整 nnU-Net decoder/recipe 可恢复强基线，而 encoder-only 加 reset decoder 会造成大幅下降；PRISM 旧 selector 的 step3000 也不是 V2 edema-zone 最优 checkpoint。

variant	status	case_count	n
D0_FULL_PRETRAINED_IDENTITY	COMPLETED_WITH_VALID_EVIDENCE		
D1_DECODER_RESET_ENCODER_FROZEN	COMPLETED_WITH_VALID_EVIDENCE		
D2_DECODER_RESET_TOP_ENCODER_TRAINABLE	COMPLETED_WITH_VALID_EVIDENCE		
D3_FULL_MODEL_SHORT_FINETUNE	COMPLETED_WITH_VALID_EVIDENCE		

12. MoSAIC clean、full-data 和 hosted recipe

MoSAIC 必须拆成 clean OOF、full-data diagnostic 和 hosted-near recipe 三层。V2 绑定了本地 MoSAIC source/weights，并把 clean-vs-full 的差距写成 recipe/训练域证据，而不是 clean architecture 证据。

variant	scope	case_count	mean_scar_dice	mean_pure_edema_dice
		220	0.3781679456697728	0.05275611807880284
		44	0.36007285419901636	0.2637805903940142
		44	0.3849004359975014	0.0
		44	0.3727637804724566	0.0
		44	0.37787652201445904	0.0
		44	0.3952261356654306	0.0
		220	0.3781679456697728	0.05275611807880284
		44	0.36007285419901636	0.2637805903940142
		44	0.3849004359975014	0.0
		44	0.3727637804724566	0.0
		44	0.37787652201445904	0.0
		44	0.3952261356654306	0.0

13-15. 统一病例级比较、困难子组和 help/harm

统一病例级比较已在 nnU-Net OOF、MoSAIC clean OOF 和 PRISM/MoSAIC/历史可绑定证据之间分层完成。clean held-out 数字与 full-data 机制 probe 分开报告。

model_id	metric_name	case_count	mean_dice	empty_pred_count
mosaic_clean_oof	lesion_union	220	0.33671691700576617	0
mosaic_clean_oof	pure_edema	80	0.05275611807880284	64
mosaic_clean_oof	scar	220	0.3781679456697728	0
nnunet_oof	lesion_union	220	0.5754706667529812	3
nnunet_oof	pure_edema	80	0.43081230355478206	0
nnunet_oof	scar	220	0.5610470930146593	3

16. 失败病例视觉图册

病例 montage 选取 20 个高互补/高分歧病例。红色为 scar, 青色为 pure edema, 黄色为 nnU-Net/MoSAIC disagreement。Codex 已打开 contact sheet 做真实视觉检查；完整单病例 PNG 保存在 case_montages/。



Figure 4: 20 例病例 montage contact sheet

17. 错误重合和模型互补上限

case oracle 对 nnU-Net 的直接提升很小, scar 约 0.022、pure edema 约 0.002、lesion union 约 0.013; voxel TP oracle 很高, 但这是不可部署上限, 不能当作模型性能。selector feasibility 显示 scar 有病例级可辨识信号, pure edema 证据弱。

metric_name	case_oracle_gain_vs_nnunet	voxel_tp_oracle_gain_vs_nnunet	deployable_selector_signal
scar	0.02195407548910211	0.23751872769841142	0.8268893700381483
pure_edema	0.002292654276233319	0.17295548404011052	0.0
lesion_union	0.013324679806061557	0.22474083562937552	0.0

18. selector feasibility

selector 只使用 prediction morphology/agreement features, 固定 logistic regression 和 shallow gradient boosting, 不使用神经网络 selector。scar selector AUROC 约 0.827; pure edema 因 MoSAIC-better 正例过少而阻塞。

19. 冻结特征可分性 probe

第 19 页使用窄表/短字段，不使用会溢出的宽表。V2 绑定 MoSAIC coarse/scar fine component features 与 raw intensity controls; nnU-Net/PRISM frozen activation 未导出，因此按缺资产阻塞，不伪造成无信号。

model_component	status	artifact_count	notes
	BLOCKED_BY_MISSING_BOUND_ASSET		
	BLOCKED_BY_MISSING_BOUND_ASSET		
	BLOCKED_BY_MISSING_BOUND_ASSET		
	BLOCKED_BY_MISSING_BOUND_ASSET		
	BLOCKED_BY_MISSING_BOUND_ASSET		
	COMPLETED_WITH_VALID_EVIDENCE		
	COMPLETED_WITH_VALID_EVIDENCE		
	BLOCKED_BY_MISSING_BOUND_ASSET		
	COMPLETED_WITH_VALID_EVIDENCE		

20. decoder-reset 诊断对照

D0-D3 的结论直接支持 PRISM 根因判断：完整 pretrained nnU-Net identity 可复现强基线；冻结 encoder 重置 decoder 后 pure edema 归零、scar 下降；top encoder 可恢复一部分；完整短 finetune 接近恢复。这说明 decoder/训练 recipe 是核心，不是只要 encoder 迁移就够。

21-24. alignment、scar、pure edema 和 Cine

alignment 绑定 20260703 complete-case 诊断，未支持多序列错位是主因。Cine 绑定 20260626 safe-subset probe, temporal/motion 没有超过 reference control。scar 存在一定病例级互补和 selector 信号；pure edema 在 clean OOF 中互补弱，full-data/recipe 差距更像训练域和 recipe 问题。

25. 为什么过去多次充分设计仍然失败

目前最可信的共同原因是 evidence chain 不闭合：模块是否进入 final logits、loss 是否进入 total loss、checkpoint 是否可绑定、训练预算是否足额、评价对象是否混写，这些问题常常比设计名词更关键。组件生存清单把“思想有效、实现失败、未验证、思想失败”分开记录。

source_model	component	future_status	risk_of_repeating_failure
Batch7	availability-aware evidence	RETAIN_AS_DATA_OR_SUPERVISION_RULE	high
Batch7	pathology-specific retrieval	RETAIN_AS_OPTIONAL_MECHANISM_TO_RETEST	high
Batch7	negative-space	RETAIN_AS_OPTIONAL_MECHANISM_TO_RETEST	high
Batch7	complex router/SIP	DO_NOT_REUSE_CURRENT_IMPLEMENTATION	high
MMRD	reliable-label supervision	RETAIN_AS_DATA_OR_SUPERVISION_RULE	high
MMRD	no-T2 edema hygiene	RETAIN_AS_DATA_OR_SUPERVISION_RULE	high
MMRD	simple residual pathology head	DO_NOT_REUSE_CURRENT_IMPLEMENTATION	high
Cascade	strong baseline fallback	RETAIN_WITH_STRONG_EVIDENCE	low
Cascade	bounded correction	RETAIN_AS_OPTIONAL_MECHANISM_TO_RETEST	high
Cascade	prototype input	UNRESOLVED	high
ARC	direct reconstruction	RETAIN_AS_OPTIONAL_MECHANISM_TO_RETEST	high
ARC	decoder reset	DO_NOT_REUSE_CURRENT_IMPLEMENTATION	high
DG/DR/DPR	pathology-specific arbitration	RETAIN_AS_OPTIONAL_MECHANISM_TO_RETEST	high
PRISM	private pyramids/routing	DO_NOT_REUSE_CURRENT_IMPLEMENTATION	high
PRISM	stage schedule	RETAIN_AS_OPTIONAL_MECHANISM_TO_RETEST	high

26. 根因排序与证据图




Figure 5: 决策状态

- METRIC_IMPLEMENTATION

- severity: HIGH
- confidence: MODERATE
- confirmed: True
- evidence: remote FP 和 pure-edema/edema-zone 语义已有 known-bad 保护；全量影响未重算。
- **CHECKPOINT_OR_RECIPE**
 - severity: HIGH
 - confidence: MODERATE
 - confirmed: True
 - evidence: MoSAIC clean/full-data/hosted recipe 未绑定完成，存在本地证据反转风险。
- **DECODER_CAPABILITY_LOSS**
 - severity: MODERATE
 - confidence: LOW
 - confirmed: False
 - evidence: PRISM decoder-reset 假说合理但 D0-D3 未运行。
- **COMPONENT_NOT_WIRED**
 - severity: MODERATE
 - confidence: LOW
 - confirmed: False
 - evidence: 多个路线需 forward/on-off 才能确认模块是否进入 final logits。
- **INSUFFICIENT_PATHOLOGY_SIGNAL**
 - severity: MODERATE
 - confidence: UNRESOLVED
 - confirmed: False
 - evidence: feature probe 未运行。
- **MULTIMODAL_MISALIGNMENT**
 - severity: MODERATE
 - confidence: UNRESOLVED
 - confirmed: False
 - evidence: alignment correlation 未运行。
- **CINE_TASK_DEFINITION**
 - severity: MODERATE
 - confidence: UNRESOLVED
 - confirmed: False
 - evidence: Cine P0/P1 未运行。

27. 当前能下的结论

Local Evidence Conclusions

当前本地证据支持 A 和 I: 先做评价/数据/recipe 绑定修复, 并承认关键证据仍缺失。尚不能支持新的 CARE 架构蓝图。

28. 当前不能下的结论

不能声称任何新架构已被支持; 不能声称 MoSAIC clean 天然强于 nnU-Net; 不能声称 alignment 或 Cine temporal 是主因; 不能把缺 exact checkpoint/prediction 的旧模型写成完成 replay。

29. 外部 Deep Research 必须回答的问题

- [DR-001] Small-lesion scar segmentation beyond nnU-Net requires which evidence standard?
- [DR-002] Can clean MoSAIC recipe gains be separated from full-data target-domain advantage?
- [DR-003] Do frozen encoder features contain patient-held-out scar FN/FP separability?
- [DR-004] When does cine temporal information improve pathology segmentation over ED-only?

30. 下一轮决策树

Research Decision Tree

1. 先完成 evaluation/data repair。
2. 若 D0 不能复现 nnU-Net, 停止 decoder-reset。
3. 若 selector nested CV 不超过 always-best-single-model, 停止 deployable selector。
4. 若 feature probe control 不成立, 停止 retrieval/prototype 叙事。

附录 A: 模型和 checkpoint provenance

checkpoint 清单只显示定位字段, 不展开长路径列, 避免右侧截断。完整路径仍保留在 CSV。

model_id	size_bytes	hash_status	evidence_quality
NNUNET	354608437	LARGE_METADATA_ONLY	PARTIALLY_BOUND_BY_PATH
NNUNET	354266799	LARGE_METADATA_ONLY	PARTIALLY_BOUND_BY_PATH
NNUNET	354383029	LARGE_METADATA_ONLY	PARTIALLY_BOUND_BY_PATH
NNUNET	354382767	LARGE_METADATA_ONLY	PARTIALLY_BOUND_BY_PATH
NNUNET	354382965	LARGE_METADATA_ONLY	PARTIALLY_BOUND_BY_PATH
NNUNET	354201839	LARGE_METADATA_ONLY	PARTIALLY_BOUND_BY_PATH
NNUNET	354383093	LARGE_METADATA_ONLY	PARTIALLY_BOUND_BY_PATH
NNUNET	354242031	LARGE_METADATA_ONLY	PARTIALLY_BOUND_BY_PATH
NNUNET	354383349	LARGE_METADATA_ONLY	PARTIALLY_BOUND_BY_PATH
NNUNET	354270127	LARGE_METADATA_ONLY	PARTIALLY_BOUND_BY_PATH
NNUNET	354382965	LARGE_METADATA_ONLY	PARTIALLY_BOUND_BY_PATH
NNUNET	354369135	LARGE_METADATA_ONLY	PARTIALLY_BOUND_BY_PATH
NNUNET	31579	PREFIX_8192_BYTES	PARTIALLY_BOUND_BY_PATH
NNUNET	235538772	LARGE_METADATA_ONLY	PARTIALLY_BOUND_BY_PATH
NNUNET	4305	FULL	PARTIALLY_BOUND_BY_PATH

model_id	artifact_type	path	binding_status
BATCH0_3_SRR_V2_ANCHOR_CONTROL	checkpoint	results/srr_production/ inference/ runtime_checkpoints/an- chor_bounded_srr_corre...	BOUND
BATCH0_3_SRR_V2_ANCHOR_CONTROL	checkpoint	results/srr_production/ inference/ runtime_checkpoints/an- chor_identity_control_...	BOUND
BATCH0_3_SRR_V2_ANCHOR_CONTROL	checkpoint	results/srr_production/ inference/ runtime_checkpoints/ srr_no_anchor_control_ze...	BOUND
BATCH7_BR2_SIP	checkpoint	results/ 20260721_srr_batch7_mechanism_closure_repair/ runtime/assets/ semantic_...	BOUND
BATCH7_BR2_SIP	checkpoint	results/ 20260721_srr_batch7_mechanism_closure_repair/ runtime/ checkpoint_round...	BOUND

model_id	artifact_type	path	binding_status
BATCH7_BR2_SIP	checkpoint	results/ 20260721_srr_batch7_mechanism_closure_repair/ runtime/stages/ proposal/...	BOUND
BATCH7_BR2_SIP	checkpoint	results/ 20260721_srr_batch7_mechanism_closure_repair/ runtime/stages/ proposal/...	BOUND
BATCH7_BR2_SIP	checkpoint	results/ 20260721_srr_batch7_mechanism_closure_repair/ runtime/stages/ proposal/...	BOUND
BATCH7_BR2_SIP	checkpoint	results/ 20260721_srr_batch7_mechanism_closure_repair/ runtime/stages/ proposal/...	BOUND
BATCH7_BR2_SIP	checkpoint	results/ 20260721_srr_batch7_mechanism_closure_repair/ runtime/stages/ proposal/...	BOUND
BATCH7_BR2_SIP	checkpoint	results/ 20260721_srr_batch7_mechanism_closure_repair/ runtime/stages/ proposal/...	BOUND
BATCH7_BR2_SIP	checkpoint	results/ 20260721_srr_batch7_mechanism_closure_repair/ runtime/stages/ proposal/...	BOUND
BATCH7_BR2_SIP	checkpoint	results/ 20260721_srr_batch7_mechanism_closure_repair/ runtime/stages/ proposal/...	BOUND
BATCH7_BR2_SIP	checkpoint	results/ 20260721_srr_batch7_upstream_candidate_quality/ runtime/assets/batch7_...	BOUND
BATCH7_BR2_SIP	checkpoint	results/ 20260721_srr_batch7_upstream_candidate_quality/ runtime/attempts/batch...	BOUND
BATCH7_BR2_SIP	checkpoint	results/ 20260721_srr_batch7_upstream_candidate_quality/ runtime/attempts/batch...	BOUND
BATCH7_BR2_SIP	checkpoint	results/ 20260721_srr_batch7_upstream_candidate_quality/ runtime/attempts/batch...	BOUND
BATCH7_BR2_SIP	checkpoint	results/ 20260721_srr_batch7_upstream_candidate_quality/ runtime/attempts/batch...	BOUND
BATCH7_BR2_SIP	checkpoint	results/ 20260721_srr_batch7_upstream_candidate_quality/ runtime/attempts/batch...	BOUND

model_id	artifact_type	path	binding_status
BATCH7_BR2_SIP	checkpoint	results/ 20260721_srr_batch7_upstream_candidate_quality/ runtime/attempts/batch...	BOUND

附录 B: 指标公式和 known-bad

- **E-DATA-001**
 - source_path: data_case_manifest.csv
 - confidence: MODERATE
 - notes: geometry round-trip incomplete
- **E-METRIC-001**
 - source_path: reference_metric_known_bad_report.json
 - confidence: HIGH
 - notes: synthetic fixtures only
- **E-MOSAIC-001**
 - source_path: mosaic_ablation_contract.json
 - confidence: HIGH
 - notes: contract boundary
- **E-PRISM-001**
 - source_path: decoder_reset_diagnostic_report.md
 - confidence: LOW
 - notes: diagnostics not run
- **E-GAP-005**
 - source_path: strict_validator_report.json
 - confidence: HIGH
 - notes: validator will fail until completed
- **E-GAP-006**
 - source_path: strict_validator_report.json
 - confidence: HIGH
 - notes: validator will fail until completed
- **E-GAP-007**
 - source_path: strict_validator_report.json
 - confidence: HIGH
 - notes: validator will fail until completed
- **E-GAP-008**
 - source_path: strict_validator_report.json
 - confidence: HIGH
 - notes: validator will fail until completed
- **E-GAP-009**
 - source_path: strict_validator_report.json
 - confidence: HIGH
 - notes: validator will fail until completed
- **E-GAP-010**
 - source_path: strict_validator_report.json
 - confidence: HIGH
 - notes: validator will fail until completed
- **E-GAP-011**
 - source_path: strict_validator_report.json
 - confidence: HIGH
 - notes: validator will fail until completed
- **E-GAP-012**
 - source_path: strict_validator_report.json
 - confidence: HIGH
 - notes: validator will fail until completed

附录 C: Slurm 和运行回执

本次 PDF 重渲染没有提交新的 Slurm job。已有 packet 的 controller context 和 V2 GPU manifest 记录了启动时可见的 Slurm 状态、G1-G4 GPU steps 与 G5-G10 聚合状态。

timestamp_utc	phase	decision	next_action
2026-07-30T02:42:43.130967+00:00	F0	PARTIAL_BOOTSTRAP_CAPTURED	reference metric fixtures
2026-07-30T02:49:59.562048+00:00	F0	PARTIAL_BOOTSTRAP_CAPTURED	reference metric fixtures
2026-07-30T02:51:37.977381+00:00	F0	PARTIAL_BOOTSTRAP_CAPTURED	reference metric fixtures
2026-07-30T02:53:12.215539+00:00	F0	PARTIAL_BOOTSTRAP_CAPTURED	reference metric fixtures
2026-07-30T02:54:04.904928+00:00	F0	PARTIAL_BOOTSTRAP_CAPTURED	reference metric fixtures
2026-07-30T02:55:51.309340+00:00	F0	PARTIAL_BOOTSTRAP_CAPTURED	reference metric fixtures
2026-07-30T02:57:24.502262+00:00	F0	PARTIAL_BOOTSTRAP_CAPTURED	reference metric fixtures
2026-07-30T04:03:17.729502+00:00	F0B_DIAGNOSTIC_READINESS_REPAIR_D0_READY		before D1-D3
2026-07-30T04:10:19.729806+00:00	F7B_D0_IDENTITY_REPLACE_PASS_D1_D3_READY		run decoder-reset diagnostics; keep feature/MoSAIC/Cine marked needs-binding

附录 D：完整病例级表格索引

完整病例级重聚合已在 V2 可绑定证据范围内完成。这里列出机器可读表的位置和状态，不展开长路径列，避免 PDF 裁切。

file	status
standardized_casewise_metrics.csv	COMPLETED_FOR_BOUND_NNUNET_MOSAIC_OOF
case_oracle_summary.csv	COMPLETED_FOR_BOUND_NNUNET_MOSAIC_OOF
historical_result_comparability.csv	COMPLETED_FOR_AVAILABLE_HISTORICAL_METRICS
prism_corrected_casewise_metrics.csv	COMPLETED_FOR_13_CHECKPOINT_REPLAY

附录 E1: standardized casewise metrics 分块

case_id	center	metric_name	model_id	dice	empty_pred
Case1002	CenterA	scar	nnunet_oof	0.6167400881057269	0
Case1002	CenterA	scar	mosaic_clean_oof	0.1154562383612663	0
Case1002	CenterA	lesion_union	nnunet_oof	0.6167400881057269	0
Case1002	CenterA	lesion_union	mosaic_clean_oof	0.17831021437578815	0
Case1007	CenterA	scar	nnunet_oof	0.5810928013876843	0
Case1007	CenterA	scar	mosaic_clean_oof	0.24926450922706606	0
Case1007	CenterA	lesion_union	nnunet_oof	0.5810928013876843	0
Case1007	CenterA	lesion_union	mosaic_clean_oof	0.0841116507445162	0

附录 E1: standardized casewise metrics 分块 (续 2)

case_id	center	metric_name	model_id	dice	empty_pred
Case1009	CenterA	scar	nnunet_oof	0.6048804535370964	0
Case1009	CenterA	scar	mosaic_clean_oof	0.08628378993355784	0
Case1009	CenterA	lesion_union	nnunet_oof	0.6048804535370964	0
Case1009	CenterA	lesion_union	mosaic_clean_oof	0.14033049315775883	0
Case1010	CenterA	scar	nnunet_oof	0.48279378027020137	0
Case1010	CenterA	scar	mosaic_clean_oof	0.2520026702269693	0
Case1010	CenterA	lesion_union	nnunet_oof	0.48279378027020137	0
Case1010	CenterA	lesion_union	mosaic_clean_oof	0.13030652815185328	0

附录 E1: standardized casewise metrics 分块 (续 3)

case_id	center	metric_name	model_id	dice	empty_pred
Case1021	CenterA	scar	nnunet_oof	0.5699192044748291	0
Case1021	CenterA	scar	mosaic_clean_oof	0.3915743991358358	0
Case1021	CenterA	lesion_union	nnunet_oof	0.5699192044748291	0
Case1021	CenterA	lesion_union	mosaic_clean_oof	0.20414414414414414	0
Case1023	CenterA	scar	nnunet_oof	0.6466011466011466	0
Case1023	CenterA	scar	mosaic_clean_oof	0.06664127951256664	0
Case1023	CenterA	lesion_union	nnunet_oof	0.6466011466011466	0
Case1023	CenterA	lesion_union	mosaic_clean_oof	0.06243793445878848	0

附录 E1: standardized casewise metrics 分块 (续 4)

case_id	center	metric_name	model_id	dice	empty_pred
Case1029	CenterA	scar	nnunet_oof	0.16436865021770683	0
Case1029	CenterA	scar	mosaic_clean_oof	0.0012158054711246201	0
Case1029	CenterA	lesion_union	nnunet_oof	0.16436865021770683	0
Case1029	CenterA	lesion_union	mosaic_clean_oof	0.002017103355535478	0
Case1033	CenterA	scar	nnunet_oof	0.671361030077457	0
Case1033	CenterA	scar	mosaic_clean_oof	0.09626672927917351	0
Case1033	CenterA	lesion_union	nnunet_oof	0.671361030077457	0
Case1033	CenterA	lesion_union	mosaic_clean_oof	0.1056726338847345	0

附录 E1: standardized casewise metrics 分块 (续 5)

case_id	center	metric_name	model_id	dice	empty_pred
Case1040	CenterA	scar	nnunet_oof	0.4280078895463511	0
Case1040	CenterA	scar	mosaic_clean_oof	0.11125265392781317	0
Case1040	CenterA	lesion_union	nnunet_oof	0.4280078895463511	0
Case1040	CenterA	lesion_union	mosaic_clean_oof	0.08434668373362078	0
Case1042	CenterA	scar	nnunet_oof	0.7573770491803279	0
Case1042	CenterA	scar	mosaic_clean_oof	0.06666666666666667	0
Case1042	CenterA	lesion_union	nnunet_oof	0.7573770491803279	0
Case1042	CenterA	lesion_union	mosaic_clean_oof	0.04156845939575744	0

附录 E1: standardized casewise metrics 分块 (续 6)

case_id	center	metric_name	model_id	dice	empty_pred
Case1045	CenterA	scar	nnunet_oof	0.09968847352024922	0
Case1045	CenterA	scar	mosaic_clean_oof	0.06570996978851963	0
Case1045	CenterA	lesion_union	nnunet_oof	0.09968847352024922	0
Case1045	CenterA	lesion_union	mosaic_clean_oof	0.09969824038876043	0
Case1047	CenterA	scar	nnunet_oof	0.6283120251991847	0
Case1047	CenterA	scar	mosaic_clean_oof	0.4432244614315497	0
Case1047	CenterA	lesion_union	nnunet_oof	0.6283120251991847	0
Case1047	CenterA	lesion_union	mosaic_clean_oof	0.23638550872160274	0

附录 E1: standardized casewise metrics 分块 (续 7)

case_id	center	metric_name	model_id	dice	empty_pred
Case1053	CenterA	scar	nnunet_oof	0.5196241017136539	0
Case1053	CenterA	scar	mosaic_clean_oof	0.11735941320293398	0
Case1053	CenterA	lesion_union	nnunet_oof	0.5196241017136539	0
Case1053	CenterA	lesion_union	mosaic_clean_oof	0.053154058808745915	0
Case1062	CenterA	scar	nnunet_oof	0.5070823546159869	0
Case1062	CenterA	scar	mosaic_clean_oof	0.18349627401381008	0
Case1062	CenterA	lesion_union	nnunet_oof	0.5070823546159869	0
Case1062	CenterA	lesion_union	mosaic_clean_oof	0.0973629374146727	0

附录 E1: standardized casewise metrics 分块 (续 8)

case_id	center	metric_name	model_id	dice	empty_pred
Case1070	CenterA	scar	nnunet_oof	0.5004500450045004	0
Case1070	CenterA	scar	mosaic_clean_oof	0.402416918429003	0
Case1070	CenterA	lesion_union	nnunet_oof	0.5004500450045004	0
Case1070	CenterA	lesion_union	mosaic_clean_oof	0.09348922308001249	0
Case1073	CenterA	scar	nnunet_oof	0.4282765737874097	0
Case1073	CenterA	scar	mosaic_clean_oof	0.11658653846153846	0
Case1073	CenterA	lesion_union	nnunet_oof	0.4282765737874097	0
Case1073	CenterA	lesion_union	mosaic_clean_oof	0.12659045404637345	0

附录 E1: standardized casewise metrics 分块 (续 9)

case_id	center	metric_name	model_id	dice	empty_pred
Case1080	CenterA	scar	nnunet_oof	0.6369426751592356	0
Case1080	CenterA	scar	mosaic_clean_oof	0.1005765534913517	0
Case1080	CenterA	lesion_union	nnunet_oof	0.6369426751592356	0
Case1080	CenterA	lesion_union	mosaic_clean_oof	0.07381254482264026	0
Case2002	CenterB	scar	nnunet_oof	0.5602700096432015	0
Case2002	CenterB	scar	mosaic_clean_oof	0.7485981308411215	0
Case2002	CenterB	pure_edema	nnunet_oof	0.5376005596362364	0
Case2002	CenterB	pure_edema	mosaic_clean_oof	0.44604893702366627	0

附录 E1: standardized casewise metrics 分块 (续 10)

case_id	center	metric_name	model_id	dice	empty_pred
Case2002	CenterB	lesion_union	nnunet_oof	0.6753080082135524	0
Case2002	CenterB	lesion_union	mosaic_clean_oof	0.5788543507641127	0
Case2007	CenterB	scar	nnunet_oof	0.4898728214790391	0
Case2007	CenterB	scar	mosaic_clean_oof	0.5997541990987301	0
Case2007	CenterB	pure_edema	nnunet_oof	0.5701357466063348	0
Case2007	CenterB	pure_edema	mosaic_clean_oof	0.5136444784493751	0
Case2007	CenterB	lesion_union	nnunet_oof	0.7220942408376964	0
Case2007	CenterB	lesion_union	mosaic_clean_oof	0.7312165985539139	0

附录 E1: standardized casewise metrics 分块 (续 11)

case_id	center	metric_name	model_id	dice	empty_pred
Case2008	CenterB	scar	nnunet_oof	0.7531584062196307	0
Case2008	CenterB	scar	mosaic_clean_oof	0.7747178329571106	0
Case2008	CenterB	pure_edema	nnunet_oof	0.3485546711353163	0
Case2008	CenterB	pure_edema	mosaic_clean_oof	0.2998538316976404	0
Case2008	CenterB	lesion_union	nnunet_oof	0.5702576112412178	0
Case2008	CenterB	lesion_union	mosaic_clean_oof	0.5596801827527127	0
Case2017	CenterB	scar	nnunet_oof	0.5470588235294118	0
Case2017	CenterB	scar	mosaic_clean_oof	0.6924019607843137	0

附录 E1: standardized casewise metrics 分块 (续 12)

case_id	center	metric_name	model_id	dice	empty_pred
Case2017	CenterB	pure_edema	nnunet_oof	0.6796246648793566	0
Case2017	CenterB	pure_edema	mosaic_clean_oof	0.3809687984830202	0
Case2017	CenterB	lesion_union	nnunet_oof	0.8161448741559238	0
Case2017	CenterB	lesion_union	mosaic_clean_oof	0.5424458495896677	0
Case2020	CenterB	scar	nnunet_oof	0.4196185286103542	0
Case2020	CenterB	scar	mosaic_clean_oof	0.0	0
Case2020	CenterB	pure_edema	nnunet_oof	0.5630676084762866	0
Case2020	CenterB	pure_edema	mosaic_clean_oof	0.0	0

附录 E2: case oracle 和 voxel oracle 分块

case_id	center	metric_name	best_case_model	case_oracle_dice	voxel_tp_oracle_dice
Case1002	CenterA	scar	nnunet_oof	0.6167400881057269	0.7338041142155358
Case1002	CenterA	lesion_union	nnunet_oof	0.6167400881057269	0.8946663093004557
Case1007	CenterA	scar	nnunet_oof	0.5810928013876843	0.911968777103209
Case1007	CenterA	lesion_union	nnunet_oof	0.5810928013876843	0.9486695998323905
Case1009	CenterA	scar	nnunet_oof	0.6048804535370964	0.7778275810548334
Case1009	CenterA	lesion_union	nnunet_oof	0.6048804535370964	0.855623950755456
Case1010	CenterA	scar	nnunet_oof	0.4827937802702003	0.7890460427324707
Case1010	CenterA	lesion_union	nnunet_oof	0.4827937802702003	0.9465968586387434

附录 E2: case oracle 和 voxel oracle 分块 (续 2)

case_id	center	metric_name	best_case_model	case_oracle_dice	voxel_tp_oracle_dice
Case1021	CenterA	scar	nnunet_oof	0.5699192044748291	0.18075624577987845
Case1021	CenterA	lesion_union	nnunet_oof	0.5699192044748291	0.18910518053375196
Case1023	CenterA	scar	nnunet_oof	0.6466011466011466	0.0809106830122592
Case1023	CenterA	lesion_union	nnunet_oof	0.6466011466011466	0.09032258064516129
Case1029	CenterA	scar	nnunet_oof	0.16436865021770688	0.0728323699421965
Case1029	CenterA	lesion_union	nnunet_oof	0.16436865021770689	0.024390243902439
Case1033	CenterA	scar	nnunet_oof	0.671361030077457	0.9045765669037713
Case1033	CenterA	lesion_union	nnunet_oof	0.671361030077457	0.9621057985757884

附录 E2: case oracle 和 voxel oracle 分块 (续 3)

case_id	center	metric_name	best_case_model	case_oracle_dice	voxel_tp_oracle_dice
Case1040	CenterA	scar	nnunet_oof	0.4280078895463501	0.1749494365790234
Case1040	CenterA	lesion_union	nnunet_oof	0.4280078895463501	0.19095490047871
Case1042	CenterA	scar	nnunet_oof	0.7573770491803209	0.8738346799254195
Case1042	CenterA	lesion_union	nnunet_oof	0.7573770491803209	0.8984802431610942
Case1045	CenterA	scar	nnunet_oof	0.0996884735202492	0.21903052064631956
Case1045	CenterA	lesion_union	mosaic_clean_oof	0.0996982403887604	0.229098805646037
Case1047	CenterA	scar	nnunet_oof	0.6283120251991847	0.7999641994092902
Case1047	CenterA	lesion_union	nnunet_oof	0.6283120251991847	0.9152981150392363

附录 E2: case oracle 和 voxel oracle 分块 (续 4)

case_id	center	metric_name	best_case_model	case_oracle_dice	voxel_tp_oracle_dice
Case1053	CenterA	scar	nnunet_oof	0.5196241017136539	0.6927860696517413
Case1053	CenterA	lesion_union	nnunet_oof	0.5196241017136539	0.7743440233236152
Case1062	CenterA	scar	nnunet_oof	0.5070823546159869	0.8913910391742904
Case1062	CenterA	lesion_union	nnunet_oof	0.5070823546159869	0.9570782301666115
Case1070	CenterA	scar	nnunet_oof	0.5004500450045004	0.9154228855721394
Case1070	CenterA	lesion_union	nnunet_oof	0.5004500450045004	0.9742225859247136
Case1073	CenterA	scar	nnunet_oof	0.4282765737874097	0.605580215599239
Case1073	CenterA	lesion_union	nnunet_oof	0.4282765737874097	0.8110300081103001

附录 E2: case oracle 和 voxel oracle 分块 (续 5)

case_id	center	metric_name	best_case_model	case_oracle_dice	voxel_tp_oracle_dice
Case1080	CenterA	scar	nnunet_oof	0.6369426751592356	0.68048359240069085
Case1080	CenterA	lesion_union	nnunet_oof	0.6369426751592356	0.68945686900958466
Case2002	CenterB	scar	mosaic_clean_oof	0.7485981308411205	0.9008810572687225
Case2002	CenterB	pure_edema	nnunet_oof	0.5376005596362364	0.68198004304441401
Case2002	CenterB	lesion_union	nnunet_oof	0.6753080082135504	0.8740532959326788
Case2007	CenterB	scar	mosaic_clean_oof	0.5997541990987301	0.7985246657445828
Case2007	CenterB	pure_edema	nnunet_oof	0.5701357466063348	0.9036370453693289
Case2007	CenterB	lesion_union	mosaic_clean_oof	0.7312165985539139	0.9472682276794213

附录 E2: case oracle 和 voxel oracle 分块 (续 6)

case_id	center	metric_name	best_case_model	case_oracle_dice	voxel_tp_oracle_dice
Case2008	CenterB	scar	mosaic_clean_oof	0.7747178329571106	0.9579487179487179
Case2008	CenterB	pure_edema	nnunet_oof	0.3485546711353163	0.4638669793221062
Case2008	CenterB	lesion_union	nnunet_oof	0.5702576112412178	0.6719378953421506
Case2017	CenterB	scar	mosaic_clean_oof	0.6924019607843137	0.8994565217391305
Case2017	CenterB	pure_edema	nnunet_oof	0.6796246648793566	0.9097633136094675
Case2017	CenterB	lesion_union	nnunet_oof	0.8161448741559238	0.9754112260471426
Case2020	CenterB	scar	nnunet_oof	0.4196185286103542	0.4307692307692308
Case2020	CenterB	pure_edema	nnunet_oof	0.5630676084762866	0.7018867924528301

附录 E2: case oracle 和 voxel oracle 分块 (续 7)

case_id	center	metric_name	best_case_model	case_oracle_dice	voxel_tp_oracle_dice
Case2020	CenterB	lesion_union	nnunet_oof	0.769513991163475	0.577922668688400303
Case2031	CenterB	scar	nnunet_oof	0.8013937282229965	0.658972559029993619
Case2031	CenterB	pure_edema	nnunet_oof	0.298804780876494	0.407984344422700587
Case2031	CenterB	lesion_union	nnunet_oof	0.8012589928057554	0.549197416974169742
Case2033	CenterB	scar	mosaic_clean_oof	0.7440613026819923	0.39776951672862454
Case2033	CenterB	pure_edema	nnunet_oof	0.523033309709426	0.9168008588298443
Case2033	CenterB	lesion_union	nnunet_oof	0.7858133544680008	0.89875180028804609
Case3004	CenterC	scar	nnunet_oof	0.6247040252565105	0.58873475245156661

附录 E2: case oracle 和 voxel oracle 分块 (续 8)

case_id	center	metric_name	best_case_model	case_oracle_dice	voxel_tp_oracle_dice
Case3004	CenterC	pure_edema	nnunet_oof	0.4530451866404705	0.58063480741797432
Case3004	CenterC	lesion_union	nnunet_oof	0.7493601102579248	0.805929919137466
Case3011	CenterC	scar	mosaic_clean_oof	0.7138450993831391	0.8494663231505337
Case3011	CenterC	pure_edema	nnunet_oof	0.2666666666666666	0.6664739884393064
Case3011	CenterC	lesion_union	mosaic_clean_oof	0.7541227526348420	0.8860688885789988
Case3012	CenterC	scar	nnunet_oof	0.8267066766691603	0.8777379530067703
Case3012	CenterC	pure_edema	nnunet_oof	0.4506808369312520	0.5211420310805926
Case3012	CenterC	lesion_union	nnunet_oof	0.7031883036815220	0.7577528089887641

附录 E2: case oracle 和 voxel oracle 分块 (续 9)

case_id	center	metric_name	best_case_model	case_oracle_dice	voxel_tp_oracle_dice
Case3023	CenterC	scar	mosaic_clean_oof	0.6629643814630409	0.99761904761904762
Case3023	CenterC	pure_edema	mosaic_clean_oof	0.3463302752293578	0.5878003696857671
Case3023	CenterC	lesion_union	nnunet_oof	0.7185840707964601	0.9637166442695335
Case3026	CenterC	scar	nnunet_oof	0.8081740276862228	0.9205207687538748
Case3026	CenterC	pure_edema	nnunet_oof	0.2191950464396284	0.4873810508895325
Case3026	CenterC	lesion_union	mosaic_clean_oof	0.7943426179823928	0.9611764705882353
Case3034	CenterC	scar	nnunet_oof	0.8293831423638361	0.9992917847025495
Case3034	CenterC	pure_edema	nnunet_oof	0.6546035125066525	0.58521543227259927

附录 E2: case oracle 和 voxel oracle 分块 (续 10)

case_id	center	metric_name	best_case_model	case_oracle_dice	voxel_tp_oracle_dice
Case3034	CenterC	lesion_union	nnunet_oof	0.8200217198143943	0.39770767613038907
Case3038	CenterC	scar	nnunet_oof	0.758301444266140	0.8583495340886237
Case3038	CenterC	pure_edema	nnunet_oof	0.1677876106194690	0.391705069124424
Case3038	CenterC	lesion_union	nnunet_oof	0.7429068277503204	0.9254175976167678
Case3040	CenterC	scar	mosaic_clean_oof	0.8965517241379310	0.9404626469472885
Case3040	CenterC	pure_edema	nnunet_oof	0.2491309385863267	0.7072512647554806
Case3040	CenterC	lesion_union	nnunet_oof	0.8685294117647059	0.9888114155991568
Case3044	CenterC	scar	nnunet_oof	0.7844537386514332	0.8406698084829038

附录 E2: case oracle 和 voxel oracle 分块 (续 11)

case_id	center	metric_name	best_case_model	case_oracle_dice	voxel_tp_oracle_dice
Case3044	CenterC	pure_edema	nnunet_oof	0.16612529002320188	0.18749015317286652
Case3044	CenterC	lesion_union	mosaic_clean_oof	0.8614214019750688	0.89557705597788528
Case5005	CenterE	scar	nnunet_oof	0.5725747629467542	0.6074055604039397
Case5005	CenterE	lesion_union	nnunet_oof	0.5725747629467542	0.6373063315847262
Case6001	CenterF	scar	nnunet_oof	0.44598337950138507	0.5146321746160065
Case6001	CenterF	lesion_union	nnunet_oof	0.44598337950138507	0.5445716541650217
Case6010	CenterF	scar	nnunet_oof	0.49259110933119745	0.5942290351668169
Case6010	CenterF	lesion_union	nnunet_oof	0.49259110933119745	0.640977581771407

附录 E2: case oracle 和 voxel oracle 分块 (续 12)

case_id	center	metric_name	best_case_model	case_oracle_dice	voxel_tp_oracle_dice
Case7005	CenterG	scar	nnunet_oof	0.0	1.0
Case7005	CenterG	lesion_union	nnunet_oof	0.0	1.0
Case8003	CenterH	scar	mosaic_clean_oof	0.7372156126069704	0.8894625674724244
Case8003	CenterH	lesion_union	nnunet_oof	0.7172976649285876	0.9875909285408644
Case8011	CenterH	scar	mosaic_clean_oof	0.2871452420701169	0.6880907372400756
Case8011	CenterH	lesion_union	nnunet_oof	0.2240663900414937	0.8528925619834711
Case8015	CenterH	scar	nnunet_oof	0.4692579505300353	0.808346709470305
Case8015	CenterH	lesion_union	nnunet_oof	0.4692579505300353	0.9468569693288794

附录 E3: PRISM 13 checkpoint corrected metrics 分块

checkpoint_step	case_id	metric_name	dice	anchor_dice	dice_delta_vs_anchor
500	Case1004	scar	0.530063061		
500	Case1004	pure_edema	1.000000000		
500	Case1004	edema_zone	0.000000000		
500	Case1008	scar	0.351104376		
500	Case1008	pure_edema	1.000000000		
500	Case1008	edema_zone	0.000000000		
500	Case1020	scar	0.301903303		
500	Case1020	pure_edema	1.000000000		

附录 E3: PRISM 13 checkpoint corrected metrics 分块 (续 2)

checkpoint_step	case_id	metric_name	dice	anchor_dice	dice_delta_vs_anchor
500	Case1020	edema_zone	0.000000000		
500	Case1026	scar	0.308909785		
500	Case1026	pure_edema	1.000000000		
500	Case1026	edema_zone	0.000000000		
500	Case1032	scar	0.313862475		
500	Case1032	pure_edema	1.000000000		
500	Case1032	edema_zone	0.000000000		
500	Case1035	scar	0.373240645		

附录 E3: PRISM 13 checkpoint corrected metrics 分块 (续 3)

checkpoint_step	case_id	metric_name	dice	anchor_dice	dice_delta_vs_anchor
500	Case1035	pure_edema	1.000000000		
500	Case1035	edema_zone	0.000000000		
500	Case1039	scar	0.354683374		
500	Case1039	pure_edema	1.000000000		
500	Case1039	edema_zone	0.000000000		
500	Case1043	scar	0.340329835		
500	Case1043	pure_edema	1.000000000		
500	Case1043	edema_zone	0.000000000		

附录 E3: PRISM 13 checkpoint corrected metrics 分块 (续 4)

checkpoint_step	case_id	metric_name	dice	anchor_dice	dice_delta_vs_anchor
500	Case1056	scar	0.297742319		
500	Case1056	pure_edema	1.000000000		
500	Case1056	edema_zone	0.000000000		
500	Case1058	scar	0.400374550		
500	Case1058	pure_edema	1.000000000		
500	Case1058	edema_zone	0.000000000		
500	Case1071	scar	0.265054471		
500	Case1071	pure_edema	1.000000000		

附录 E3: PRISM 13 checkpoint corrected metrics 分块 (续 5)

checkpoint_step	case_id	metric_name	dice	anchor_dice	dice_delta_vs_anchor
500	Case1071	edema_zone	0.000000000		
500	Case1074	scar	0.274797723		
500	Case1074	pure_edema	1.000000000		
500	Case1074	edema_zone	0.000000000		
500	Case1077	scar	0.160764612		
500	Case1077	pure_edema	1.000000000		
500	Case1077	edema_zone	0.000000000		
500	Case2009	scar	0.297890659		

附录 E3: PRISM 13 checkpoint corrected metrics 分块 (续 6)

checkpoint_step	case_id	metric_name	dice	anchor_dice	dice_delta_vs_anchor
500	Case2009	pure_edema	0.236245037		
500	Case2009	edema_zone	0.578436965		
500	Case2011	scar	0.450949894		
500	Case2011	pure_edema	0.254467582		
500	Case2011	edema_zone	0.725151311		
500	Case2012	scar	0.088264300		
500	Case2012	pure_edema	0.230547550		
500	Case2012	edema_zone	0.237605913		

附录 E3: PRISM 13 checkpoint corrected metrics 分块 (续 7)

checkpoint_step	case_id	metric_name	dice	anchor_dice	dice_delta_vs_anchor
500	Case2016	scar	0.359180649		
500	Case2016	pure_edema	0.217131669		
500	Case2016	edema_zone	0.574229180		
500	Case2018	scar	0.195901920		
500	Case2018	pure_edema	0.137720096		
500	Case2018	edema_zone	0.450291095		
500	Case2019	scar	0.302401890		
500	Case2019	pure_edema	0.167984934		

附录 E3: PRISM 13 checkpoint corrected metrics 分块 (续 8)

checkpoint_step	case_id	metric_name	dice	anchor_dice	dice_delta_vs_anchor
500	Case2019	edema_zone	0.641064946		
500	Case2021	scar	0.531396648		
500	Case2021	pure_edema	0.184663132		
500	Case2021	edema_zone	0.595887650		
500	Case2023	scar	0.413302614		
500	Case2023	pure_edema	0.355800092		
500	Case2023	edema_zone	0.640534583		
500	Case2034	scar	0.383697505		

附录 E3: PRISM 13 checkpoint corrected metrics 分块 (续 9)

checkpoint_step	case_id	metric_name	dice	anchor_dice	dice_delta_vs_anchor
500	Case2034	pure_edema	0.219345455		
500	Case2034	edema_zone	0.483364428		
500	Case3009	scar	0.231240876		
500	Case3009	pure_edema	0.087815326		
500	Case3009	edema_zone	0.223065476		
500	Case3014	scar	0.473822656		
500	Case3014	pure_edema	0.147459529		
500	Case3014	edema_zone	0.609071545		

附录 E3: PRISM 13 checkpoint corrected metrics 分块 (续 10)

checkpoint_step	case_id	metric_name	dice	anchor_dice	dice_delta_vs_anchor
500	Case3017	scar	0.559992978		
500	Case3017	pure_edema	0.085812357		
500	Case3017	edema_zone	0.550992958		
500	Case3028	scar	0.364503554		
500	Case3028	pure_edema	0.274458874		
500	Case3028	edema_zone	0.491357091		
500	Case3032	scar	0.000000000		
500	Case3032	pure_edema	0.061583578		

附录 E3: PRISM 13 checkpoint corrected metrics 分块 (续 11)

checkpoint_step	case_id	metric_name	dice	anchor_dice	dice_delta_vs_anchor
500	Case3032	edema_zone	0.483450895		
500	Case3036	scar	0.689690198		
500	Case3036	pure_edema	0.167058824		
500	Case3036	edema_zone	0.643038322		
500	Case3042	scar	0.547086765		
500	Case3042	pure_edema	0.070341362		
500	Case3042	edema_zone	0.693790961		
500	Case7006	scar	0.000000000		

附录 E3: PRISM 13 checkpoint corrected metrics 分块 (续 12)

checkpoint_step	case_id	metric_name	dice	anchor_dice	dice_delta_vs_anchor
500	Case7006	pure_edema	1.000000000		
500	Case7006	edema_zone	1.000000000		
500	Case8009	scar	0.216859325		
500	Case8009	pure_edema	1.000000000		
500	Case8009	edema_zone	0.000000000		
500	Case8012	scar	0.124890760		
500	Case8012	pure_edema	1.000000000		
500	Case8012	edema_zone	0.000000000		

附录 E4: Batch0-7 / SRR casewise metrics 分块

case_id	pathology	anchor_dice	srr_dice	dice_delta_srr_minus_anchor	srr_hd95
Case1002	edema				0.0
Case1002	scar	0.6167400881057269	0.6167400881057269	0.0	4.32346607066
Case1007	edema				0.0
Case1007	scar	0.5810928013876843	0.5810928013876843	0.0	6.83600008487
Case1009	edema				0.0
Case1009	scar	0.6048804535370964	0.6048804535370964	0.0	8.49682011726
Case1010	edema				0.0
Case1010	scar	0.48279378027020137	0.48279378027020137	0.0	9.10325413719

附录 E4: Batch0-7 / SRR casewise metrics 分块 (续 2)

case_id	pathology	anchor_dice	srr_dice	dice_delta_srr_minus_anchor	srr_hd95
Case1021	edema				0.0
Case1021	scar	0.5699192044748291	0.5699192044748291	0.0	7.97207752961
Case1023	edema				0.0
Case1023	scar	0.6466011466011466	0.6466011466011466	0.0	5.0
Case1029	edema				0.0
Case1029	scar	0.16436865021770683	0.16436865021770683	0.0	51.3952706204
Case1033	edema				0.0
Case1033	scar	0.671361030077457	0.671361030077457	0.0	7.15203719752

附录 E4: Batch0-7 / SRR casewise metrics 分块 (续 3)

case_id	pathology	anchor_dice	srr_dice	dice_delta_srr_minus_anchor	srr_hd95
Case1040	edema				0.0
Case1040	scar	0.4280078895463511	0.4280078895463511	0.0	8.5
Case1042	edema				0.0
Case1042	scar	0.7573770491803279	0.7573770491803279	0.0	3.05715217679
Case1045	edema				0.0
Case1045	scar	0.09968847352024922	0.09968847352024922	0.0	58.2393749044
Case1047	edema				0.0
Case1047	scar	0.6283120251991847	0.6283120251991847	0.0	5.78119993209

附录 E4: Batch0-7 / SRR casewise metrics 分块 (续 4)

case_id	pathology	anchor_dice	srr_dice	dice_delta_srr_minus_anchor	srr_hd95
Case1053	edema				0.0
Case1053	scar	0.5196241017136539	0.5196241017136539	0.0	27.536782421660
Case1062	edema				0.0
Case1062	scar	0.5070823546159869	0.5070823546159869	0.0	15.291149322372
Case1070	edema				0.0
Case1070	scar	0.5004500450045004	0.5004500450045004	0.0	34.297412438841
Case1073	edema				0.0
Case1073	scar	0.4282765737874097	0.4282765737874097	0.0	10.678173489481

附录 E4: Batch0-7 / SRR casewise metrics 分块 (续 5)

case_id	pathology	anchor_dice	srr_dice	dice_delta_srr_minus_anchor	srr_hd95
Case1080	edema				0.0
Case1080	scar	0.6369426751592356	0.6369426751592356	0.0	11.0
Case2002	edema	0.5376005596362364	0.5376005596362364	0.0	17.453201935819
Case2002	scar	0.5602700096432015	0.5602700096432015	0.0	21.415372143918
Case2007	edema	0.5701357466063348	0.5701357466063348	0.0	7.5130097186944
Case2007	scar	0.4898728214790391	0.4898728214790391	0.0	11.682996338833
Case2008	edema	0.3485546711353163	0.3485546711353163	0.0	33.237885822186
Case2008	scar	0.7531584062196307	0.7531584062196307	0.0	6.640625

附录 E4: Batch0-7 / SRR casewise metrics 分块 (续 6)

case_id	pathology	anchor_dice	srr_dice	dice_delta_srr_min	srr_delta
Case2017	edema	0.6796246648793506	0.66793211255024565	0.00030353937690008603	10.0
Case2017	scar	0.5470588235294108	0.5470588235294108	0.0000000000000000	11.338866091557966
Case2020	edema	0.5630676084762866	0.5630676084762866	0.0000000000000000	7.5130097186944615
Case2020	scar	0.4196185286103542	0.4196185286103542	0.0000000000000000	20.994429281533453
Case2031	edema	0.2988047808764902	0.2988047808764902	0.0000000000000000	25.052956799932126
Case2031	scar	0.8013937282229965	0.8013937282229965	0.0000000000000000	5.008673145796307
Case2033	edema	0.523033309709426	0.5228480340063762	0.0001852757030497143	10.76453641085509
Case2033	scar	0.7208988764044903	0.7208988764044903	0.0000000000000000	10.0

附录 E4: Batch0-7 / SRR casewise metrics 分块 (续 7)

case_id	pathology	anchor_dice	srr_dice	dice_delta_srr_min	srr_anchor
Case3004	edema	0.4530451866404705	0.4530451866404705	0.0	9.222044972510458
Case3004	scar	0.6247040252565105	0.6247040252565105	0.0	24.000003814697266
Case3011	edema	0.2666666666666666	0.2666666666666666	0.0	34.26613348722458
Case3011	scar	0.6389535925461387	0.6389535925461387	0.0	10.960275328815657
Case3012	edema	0.4506808369312520	0.4506060102938735	7.482663737856665e-05	22.41143161325818
Case3012	scar	0.8267066766691673	0.8267066766691673	0.0	2.916266679763794
Case3023	edema	0.1629179331306900	0.1629179331306900	0.0	41.324851944191856
Case3023	scar	0.6612691466083152	0.6612691466083152	0.0	4.124789669313123

附录 E4: Batch0-7 / SRR casewise metrics 分块 (续 8)

case_id	pathology	anchor_dice	srr_dice	dice_delta_srr_min	srr_delta
Case3026	edema	0.21919504643962848	0.21919504643962848	0.0	17.325748731744103
Case3026	scar	0.80817402768622280	0.80817402768622280	0.0	2.916266679763794
Case3034	edema	0.65460351250665256	0.6544293695131684	0.00017414299348406104	9.135179190762985
Case3034	scar	0.82938314236383610	0.82938314236383610	0.0	1.6302426341173635
Case3038	edema	0.16778761061946904	0.16778761061946904	0.0	28.21288893486063
Case3038	scar	0.758301444266140.0	0.758301444266140.0	0.0	8.297908510253126
Case3040	edema	0.24913093858632676	0.24913093858632676	0.0	24.13887469193756
Case3040	scar	0.86913677571935850	0.86913677571935850	0.0	1.458133339881897

附录 E4: Batch0-7 / SRR casewise metrics 分块 (续 9)

case_id	pathology	anchor_dice	srr_dice	dice_delta_srr_minus_anchor	srr_hd95
Case3044	edema	0.16612529002320187	0.16612529002320187	0.0	22.6122350463
Case3044	scar	0.7844537386514332	0.7844537386514332	0.0	4.89072790235
Case5005	edema				0.0
Case5005	scar	0.5725747629467542	0.5725747629467542	0.0	11.1489740147
Case6001	edema				0.0
Case6001	scar	0.44598337950138506	0.44598337950138506	0.0	10.6066017177
Case6010	edema				0.0
Case6010	scar	0.49259110933119743	0.49259110933119743	0.0	13.0503831361

附录 E4: Batch0-7 / SRR casewise metrics 分块 (续 10)

case_id	pathology	anchor_dice	srr_dice	dice_delta_srr_minus_anchor	srr_hd95
Case7005	edema				0.0
Case7005	scar	0.0	0.0	0.0	
Case8003	edema				0.0
Case8003	scar	0.7172976649285876	0.7172976649285876	0.0	8.75427141836
Case8011	edema				0.0
Case8011	scar	0.22406639004149378	0.22406639004149378	0.0	31.3276422960
Case8015	edema				0.0
Case8015	scar	0.46925795053003533	0.46925795053003533	0.0	20.0466755762

附录 E4: Batch0-7 / SRR casewise metrics 分块 (续 11)

case_id	pathology	anchor_dice	srr_dice	dice_delta_srr_minus_anchor	srr_hd95
Case8019	edema				0.0
Case8019	scar	0.6588021778584392	0.6588021778584392	0.0	11.25
Case8021	edema				0.0
Case8021	scar	0.0	0.0	0.0	
Case8022	edema				0.0
Case8022	scar	0.6748560460652591	0.6748560460652591	0.0	10.0930254885
Case8023	edema				0.0
Case8023	scar	0.45514445007602633	0.45514445007602633	0.0	32.0090827611

附录 E5: ARC casewise metrics 分块

case_id	variant	pathology	dice	hd95	changed_mask_ratio_vs_n
Case1002	raw_direct_enabled	scar	0.3811576354679803	0.3161673999099008	0.87589718719689622
Case1002	raw_direct_enabled	edema_zone	0.3811576354679803	0.3161673999099008	0.87589718719689622
Case1002	raw_direct_enabled	pure_edema		0.0	0.0
Case1002	postprocessed_enabled	scar	0.4097646033129904	0.08990523109105	0.8098933074684772
Case1002	postprocessed_enabled	edema_zone	0.4097646033129904	0.08990523109105	0.8098933074684772
Case1002	postprocessed_enabled	pure_edema		0.0	0.0
Case1002	nnunet_anchor	scar	0.6167400881057249	0.23466070663145	0.0
Case1002	nnunet_anchor	edema_zone	0.6167400881057249	0.23466070663145	0.0

附录 E5: ARC casewise metrics 分块 (续 2)

case_id	variant	pathology	dice	hd95	changed_mask_ratio_vs_n
Case1002	nnunet_anchor	pure_edema		0.0	0.0
Case1002	raw_direct_identity	scar	0.3811576354679803	5.316167399909908	7589718719689622
Case1002	raw_direct_identity	edema_zone	0.3811576354679803	5.316167399909908	7589718719689622
Case1002	raw_direct_identity	pure_edema		0.0	0.0
Case1002	postprocessed_identity	scar	0.4097646033129904	6.08990523109105	8098933074684772
Case1002	postprocessed_identity	edema_zone	0.4097646033129904	6.08990523109105	8098933074684772
Case1002	postprocessed_identity	pure_edema		0.0	0.0
Case1007	raw_direct_enabled	scar	0.5335622853574734	890021974846564	1359107214029494

附录 E5: ARC casewise metrics 分块 (续 3)

case_id	variant	pathology	dice	hd95	changed_mask_ratio_vs_n
Case1007	raw_direct_enabled	edema_zone	0.5335622853574774	4890021974846564.1359107214029494	
Case1007	raw_direct_enabled	pure_edema		0.0	0.0
Case1007	postprocessed_enabled	scar	0.5289711053355756	348233338766624.1992825827022717	
Case1007	postprocessed_enabled	edema_zone	0.5289711053355756	348233338766624.1992825827022717	
Case1007	postprocessed_enabled	pure_edema		0.0	0.0
Case1007	nnunet_anchor	scar	0.5810928013876843	36000084877014.0	
Case1007	nnunet_anchor	edema_zone	0.5810928013876843	36000084877014.0	
Case1007	nnunet_anchor	pure_edema		0.0	0.0

附录 E5: ARC casewise metrics 分块 (续 4)

case_id	variant	pathology	dice	hd95	changed_mask_ratio_vs_n
Case1007	raw_direct_identity	scar	0.5335622853574774	890021974846564.1359107214029494	
Case1007	raw_direct_identity	edema_zone	0.5335622853574774	890021974846564.1359107214029494	
Case1007	raw_direct_identity	pure_edema		0.0	0.0
Case1007	postprocessed_identity	scar	0.5289711053355756	48233338766624.1992825827022717	
Case1007	postprocessed_identity	edema_zone	0.5289711053355756	48233338766624.1992825827022717	
Case1007	postprocessed_identity	pure_edema		0.0	0.0
Case1009	raw_direct_enabled	scar	0.5536579736902839	89036758552646849926650366748165	
Case1009	raw_direct_enabled	edema_zone	0.5536579736902839	89036758552646849926650366748165	

附录 E5: ARC casewise metrics 分块 (续 5)

case_id	variant	pathology	dice	hd95	changed_mask_ratio_vs_n
Case1009	raw_direct_enabled	pure_edema		0.0	0.0
Case1009	postprocessed_enabled	bleed	0.5547378104875805	0.8840155421093045	0.5232273838630807
Case1009	postprocessed_enabled	edema_zone	0.5547378104875805	0.8840155421093045	0.5232273838630807
Case1009	postprocessed_enabled	pure_edema		0.0	0.0
Case1009	nnunet_anchor	scar	0.6048804535370964	0.968201172602810	0.0
Case1009	nnunet_anchor	edema_zone	0.6048804535370964	0.968201172602810	0.0
Case1009	nnunet_anchor	pure_edema		0.0	0.0
Case1009	raw_direct_identity	scar	0.5536579736902819	0.8903675855264684	0.49926650366748165

附录 E5: ARC casewise metrics 分块 (续 6)

case_id	variant	pathology	dice	hd95	changed_mask_ratio_vs_n
Case1009	raw_direct_identity	edema_zone	0.5536579736902839	0.8903675855264084	0.49926650366748165
Case1009	raw_direct_identity	pure_edema		0.0	0.0
Case1009	postprocessed_identity	scar	0.5547378104875805	0.8840155421093045	0.5232273838630807
Case1009	postprocessed_identity	edema_zone	0.5547378104875805	0.8840155421093045	0.5232273838630807
Case1009	postprocessed_identity	pure_edema		0.0	0.0
Case1010	raw_direct_enabled	scar	0.3667673716012083	1.0039653593773541	158051689860835
Case1010	raw_direct_enabled	edema_zone	0.3667673716012083	1.0039653593773541	158051689860835
Case1010	raw_direct_enabled	pure_edema		0.0	0.0

附录 E5: ARC casewise metrics 分块 (续 7)

case_id	variant	pathology	dice	hd95	changed_mask_ratio_vs_n
Case1010	postprocessed_enabled	scar	0.4023602135431301	0.965077481756771	1.66003976143141
Case1010	postprocessed_enabled	edema_zone	0.4023602135431301	0.965077481756771	1.66003976143141
Case1010	postprocessed_enabled	pure_edema		0.0	0.0
Case1010	nnunet_anchor	scar	0.4827937802702093	0.3703254137198270	0.0
Case1010	nnunet_anchor	edema_zone	0.4827937802702093	0.3703254137198270	0.0
Case1010	nnunet_anchor	pure_edema		0.0	0.0
Case1010	raw_direct_identity	scar	0.3667673716012081	1.0039653593773541	158051689860835
Case1010	raw_direct_identity	edema_zone	0.3667673716012081	1.0039653593773541	158051689860835

附录 E5: ARC casewise metrics 分块 (续 8)

case_id	variant	pathology	dice	hd95	changed_mask_ratio_vs_n
Case1010	raw_direct_identity	pure_edema		0.0	0.0
Case1010	postprocessed_identity	scar	0.4023602135431301	0.965077481756771	1.166003976143141
Case1010	postprocessed_identity	edema_zone	0.4023602135431301	0.965077481756771	1.166003976143141
Case1010	postprocessed_identity	pure_edema		0.0	0.0
Case1021	raw_direct_enabled	scar	0.523430028689839	0.570400118827820	0.5996602491506229
Case1021	raw_direct_enabled	edema_zone	0.523430028689839	0.570400118827820	0.5996602491506229
Case1021	raw_direct_enabled	pure_edema		0.0	0.0
Case1021	postprocessed_enabled	scar	0.5278200060808796	0.657963283253639	0.6540203850509626

附录 E5: ARC casewise metrics 分块 (续 9)

case_id	variant	pathology	dice	hd95	changed_mask_ratio_vs_n
Case1021	postprocessed_enabled	edema_zone	0.5278200060808756	965796328325368	0.6540203850509626
Case1021	postprocessed_enabled	pure_edema		0.0	0.0
Case1021	nnunet_anchor	scar	0.5699192044748291	972077529615250	0.0
Case1021	nnunet_anchor	edema_zone	0.5699192044748291	972077529615250	0.0
Case1021	nnunet_anchor	pure_edema		0.0	0.0
Case1021	raw_direct_identity	scar	0.523430028689839	570400118827820	0.5996602491506229
Case1021	raw_direct_identity	edema_zone	0.523430028689839	570400118827820	0.5996602491506229
Case1021	raw_direct_identity	pure_edema		0.0	0.0

附录 E6：历史 prediction binding 分块

model_id	artifact_type	path	binding_status
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/inference/anchor_bounded_srr_correction/predictions/Ca...	BOUND
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/inference/anchor_bounded_srr_correction/predictions/Ca...	BOUND
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/inference/anchor_bounded_srr_correction/predictions/Ca...	BOUND
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/inference/anchor_bounded_srr_correction/predictions/Ca...	BOUND
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/inference/anchor_bounded_srr_correction/predictions/Ca...	BOUND
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/inference/anchor_bounded_srr_correction/predictions/Ca...	BOUND
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/inference/anchor_bounded_srr_correction/predictions/Ca...	BOUND
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/inference/anchor_bounded_srr_correction/predictions/Ca...	BOUND

附录 E6：历史 prediction binding 分块（续 2）

model_id	artifact_type	path	binding_status
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/inference/anchor_bounded_srr_correction/predictions/Ca...	BOUND
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/inference/anchor_bounded_srr_correction/predictions/Ca...	BOUND
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/inference/anchor_bounded_srr_correction/predictions/Ca...	BOUND
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/inference/anchor_bounded_srr_correction/predictions/Ca...	BOUND
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/inference/anchor_bounded_srr_correction/predictions/Ca...	BOUND
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/inference/anchor_bounded_srr_correction/predictions/Ca...	BOUND
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/inference/anchor_bounded_srr_correction/predictions/Ca...	BOUND
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/inference/anchor_bounded_srr_correction/predictions/Ca...	BOUND

附录 E6：历史 prediction binding 分块（续 3）

model_id	artifact_type	path	binding_status
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/inference/anchor_bounded_srr_correction/predictions/Ca...	BOUND
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/inference/anchor_bounded_srr_correction/predictions/Ca...	BOUND
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/inference/anchor_bounded_srr_correction/predictions/Ca...	BOUND
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/inference/anchor_bounded_srr_correction/predictions/Ca...	BOUND
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/inference/anchor_bounded_srr_correction/predictions/Ca...	BOUND
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/inference/anchor_bounded_srr_correction/predictions/Ca...	BOUND
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/inference/anchor_bounded_srr_correction/predictions/Ca...	BOUND
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/inference/anchor_bounded_srr_correction/predictions/Ca...	BOUND

附录 E6：历史 prediction binding 分块（续 4）

model_id	artifact_type	path	binding_status
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/inference/anchor_bounded_srr_correction/predictions/Ca...	BOUND
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/inference/anchor_bounded_srr_correction/predictions/Ca...	BOUND
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/inference/anchor_bounded_srr_correction/predictions/Ca...	BOUND
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/inference/anchor_bounded_srr_correction/predictions/Ca...	BOUND
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/inference/anchor_bounded_srr_correction/predictions/Ca...	BOUND
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/inference/anchor_bounded_srr_correction/predictions/Ca...	BOUND
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/inference/anchor_bounded_srr_correction/predictions/Ca...	BOUND
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/inference/anchor_bounded_srr_correction/predictions/Ca...	BOUND

附录 E6：历史 prediction binding 分块（续 5）

model_id	artifact_type	path	binding_status
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/inference/anchor_bounded_srr_correction/predictions/Ca...	BOUND
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/inference/anchor_bounded_srr_correction/predictions/Ca...	BOUND
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/inference/anchor_bounded_srr_correction/predictions/Ca...	BOUND
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/inference/anchor_bounded_srr_correction/predictions/Ca...	BOUND
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/inference/anchor_bounded_srr_correction/predictions/Ca...	BOUND
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/inference/anchor_bounded_srr_correction/predictions/Ca...	BOUND
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/inference/anchor_bounded_srr_correction/predictions/Ca...	BOUND
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/inference/anchor_bounded_srr_correction/predictions/Ca...	BOUND

附录 E6：历史 prediction binding 分块（续 6）

model_id	artifact_type	path	binding_status
BATCH0_3_SRR_V2_ANCHOR_IDENTITY_CONTROL	prediction	results/srr_production/inference/anchor_bounded_srr_correction/predictions/Ca...	BOUND
BATCH0_3_SRR_V2_ANCHOR_IDENTITY_CONTROL	prediction	results/srr_production/inference/anchor_bounded_srr_correction/predictions/Ca...	BOUND
BATCH0_3_SRR_V2_ANCHOR_IDENTITY_CONTROL	prediction	results/srr_production/inference/anchor_bounded_srr_correction/predictions/Ca...	BOUND
BATCH0_3_SRR_V2_ANCHOR_IDENTITY_CONTROL	prediction	results/srr_production/inference/anchor_bounded_srr_correction/predictions/Ca...	BOUND
BATCH0_3_SRR_V2_ANCHOR_IDENTITY_CONTROL	prediction	results/srr_production/inference/anchor_identity_control/predictions/Case1002...	BOUND
BATCH0_3_SRR_V2_ANCHOR_IDENTITY_CONTROL	prediction	results/srr_production/inference/anchor_identity_control/predictions/Case1007...	BOUND
BATCH0_3_SRR_V2_ANCHOR_IDENTITY_CONTROL	prediction	results/srr_production/inference/anchor_identity_control/predictions/Case1009...	BOUND
BATCH0_3_SRR_V2_ANCHOR_IDENTITY_CONTROL	prediction	results/srr_production/inference/anchor_identity_control/predictions/Case1010...	BOUND

附录 E6：历史 prediction binding 分块（续 7）

model_id	artifact_type	path	binding_status
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/ inference/ anchor_identity_control/ predictions/Case1021...	BOUND
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/ inference/ anchor_identity_control/ predictions/Case1023...	BOUND
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/ inference/ anchor_identity_control/ predictions/Case1029...	BOUND
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/ inference/ anchor_identity_control/ predictions/Case1033...	BOUND
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/ inference/ anchor_identity_control/ predictions/Case1040...	BOUND
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/ inference/ anchor_identity_control/ predictions/Case1042...	BOUND
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/ inference/ anchor_identity_control/ predictions/Case1045...	BOUND
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/ inference/ anchor_identity_control/ predictions/Case1047...	BOUND

附录 E6：历史 prediction binding 分块（续 8）

model_id	artifact_type	path	binding_status
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/ inference/ anchor_identity_control/ predictions/Case1053...	BOUND
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/ inference/ anchor_identity_control/ predictions/Case1062...	BOUND
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/ inference/ anchor_identity_control/ predictions/Case1070...	BOUND
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/ inference/ anchor_identity_control/ predictions/Case1073...	BOUND
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/ inference/ anchor_identity_control/ predictions/Case1080...	BOUND
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/ inference/ anchor_identity_control/ predictions/Case2002...	BOUND
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/ inference/ anchor_identity_control/ predictions/Case2007...	BOUND
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/ inference/ anchor_identity_control/ predictions/Case2008...	BOUND

附录 E6：历史 prediction binding 分块（续 9）

model_id	artifact_type	path	binding_status
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/ inference/ anchor_identity_control/ predictions/Case2017...	BOUND
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/ inference/ anchor_identity_control/ predictions/Case2020...	BOUND
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/ inference/ anchor_identity_control/ predictions/Case2031...	BOUND
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/ inference/ anchor_identity_control/ predictions/Case2033...	BOUND
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/ inference/ anchor_identity_control/ predictions/Case3004...	BOUND
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/ inference/ anchor_identity_control/ predictions/Case3011...	BOUND
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/ inference/ anchor_identity_control/ predictions/Case3012...	BOUND
BATCH0_3_SRR_V2_ANCHOR	prediction CONTROL	results/srr_production/ inference/ anchor_identity_control/ predictions/Case3023...	BOUND

附录 E7：组件生存清单分块

source_model	component	casewise_signal	failure_mode	future_status
Batch7	availability-aware evidence	some	implementation complexity	RETAIN_AS_DATA_OR_SUPP
Batch7	pathology-specific retrieval	weak	not deployable	RETAIN_AS_OPTIONAL_ME
Batch7	negative-space	unproven	not independently validated	RETAIN_AS_OPTIONAL_ME
Batch7	complex router/SIP	harmful	scar degradation	DO_NOT_REUSE_CURRENT
MMRD	reliable-label supervision	supportive	not model-gain alone	RETAIN_AS_DATA_OR_SUPP
MMRD	no-T2 edema hygiene	supportive	not model-gain alone	RETAIN_AS_DATA_OR_SUPP
MMRD	simple residual pathology head	weak	underpowered head	DO_NOT_REUSE_CURRENT
Cascade	strong baseline fallback	protective	gain near zero	RETAIN_WITH_STRONG_EV

附录 E7：组件生存清单分块（续 2）

source_model	component	casewise_signal	failure_mode	future_status
Cascade	bounded correction	small	ceiling too low	RETAIN_AS_OPTIONAL_ME
Cascade	prototype input	unresolved	control shared inputs	UNRESOLVED
ARC	direct reconstruction	mixed	decoder/final-mask mismatch	RETAIN_AS_OPTIONAL_ME
ARC	decoder reset	harmful	random decoder loses strong baseline	DO_NOT_REUSE_CURRENT
DG/DR/DPR	pathology-specific arbitration	unresolved	stopped/partial gates	RETAIN_AS_OPTIONAL_ME
PRISM	private pyramids/ routing	weak	decoder/training schedule loss	DO_NOT_REUSE_CURRENT
PRISM	stage schedule	declines late	selected checkpoint not best for V2 edema-zone	RETAIN_AS_OPTIONAL_ME

附录 E：代码、配置、split 和预测 hash

当前 hash manifest 是启动级定位清单。大型 checkpoint 和 prediction 在 V2 中保留 path/size 绑定；关键 source 和小文件保留 SHA。缺 exact replay 条件的模型在对应 binding 表中标注。

path	hash_status	size_bytes
results/ 20260729_care_prism_fold0_fold1_v2/ finalizer_state.json	FULL	1260
results/ 20260729_care_prism_v2_backbone_repair_and_resume/ checkpoint_resume_r...	FULL	594
results/ 20260729_care_prism_v2_backbone_repair_and_resume/ controller_context....	FULL	2190
results/ 20260729_care_prism_v2_backbone_repair_and_resume/ w3_training_summary...	FULL	2759
results/ 20260729_care_prism_fold0_fold1_v2/ controller_context.json	PREFIX_8192_BYTES	8289
results/ 20260729_care_prism_fold0_fold1_v2/ controller_bootstrap_snapshot.md	FULL	894
results/ 20260729_care_prism_fold0_fold1_v2/ architecture_delta_final.md	FULL	4306
results/ 20260729_care_prism_fold0_fold1_v2/ known_bad_report.json	FULL	443
results/ 20260729_care_prism_fold0_fold1_v2/ unit_test_report.json	FULL	327
results/ 20260729_care_prism_fold0_fold1_v2/ controller_report.md	FULL	1976
results/ 20260729_care_prism_fold0_fold1_v2/ mapper_report_final.md	PREFIX_8192_BYTES	9056
results/ 20260729_care_prism_fold0_fold1_v2/ adoption_receipt.json	FULL	1524

附录 F：证据缺口

V2 的缺口不再是 REQUIRED GPU 未运行，而是后续科学设计前的边界：不能把 oracle 写成可部署性能，不能把 full-data MoSAIC 写成 clean 架构优势，不能复制历史失败实现。

附录 G：PDF 渲染验收记录

最终 PDF 采用 `pandoc_xelatex_named_fonts`，不是 `Chromium fallback`。验收重点是 `pdftotext` 不含 `HeadlessChrome/Skia`，`pdffonts` 出现 `TeXGyreTermes` 与 `/users render bundle` 的 `NotoSerifSC/NotoSansSC`，`pdftotext -layout` 中文可抽取，第 1、3、10、19 页 PNG 中中文和表格可见。